



REF Memorization



Class I Medical Device

User manual

Distribution mode

Available for direct download at http://virtualisvr.com/espace-client/ Use under licence





DESCRIPTION

MEMORIZATION software is an immersive 3D simulation based on virtual reality technology, i.e. it allows a person to be immersed in an artificial digitally created world. **Memorization** is a memory assessment and training software based on the memorization of a shopping list.

INDICATIONS

Assessment and cognitive rehabilitation in the case of neurological disorders.

CONTRAINDICATIONS

Epileptic patients, children under 15 years of age, pregnant women.

FOR USE BY

Healthcare professionals: Physiotherapists; Occupational therapists; Neuropsychologists; ENT doctors; Neurologists; PMR doctors (physical medicine and rehabilitation), etc.

Research Centers: CNRS, CHU, INSERM, etc.

WARNINGS AND CAUTIONS

During sessions, stay close to the patient in order to anticipate any loss of balance or discomfort caused by the use of virtual reality.

Define a working area of about 3m² to allow for risk-free movements.

Take a 10 to 15 minute break every 30 minutes of use.

Potential adverse effects are those due to software, i.e. vomiting, malaise, dizziness, syncope

The accessories required to use the software may emit radio waves that can interfere with the operation of nearby electronic devices. If you have a pacemaker or other implanted medical device, do not use the product until you have taken advice from your doctor or the manufacturer of your medical device.



Any serious incident should be notified in writing to qualite@virtualisvr.com

VIRTUALIS Page 1 of 10



Table of Contents

1.	GEN	ERAL	3
	1.1.	Advice for use	3
	1.2.	Hardware and minimum configuration requirements	3
2.	USE	OF PATIENT MANAGEMENT	4
3.	MEN	MORIZATION	6
	3.1.	Start interface	6
	3.2.	Module field of application	7
	3.3.	Installing the patient	7
	3.4.	Session settings	7
	3.5.	Shortcuts	9
	3.6.	Data processing	. 10



1. GENERAL

1.1. Advice for use

This type of re-education must be undertaken progressively, especially in Virtual Reality where the stimulation is much more "powerful" than the traditional optokinetic stimulators.

These stimulations have the potential to cause certain disorders: Vasovagal syncope, epileptic seizures, migraines, etc. (Despite a test phase on more than 2000 patients. Similarly to previous generation optokinetics, caution is required)

The contraindications are identical: Mainly epilepsy and migraines.

As postural reactions can be spectacular, it is VERY STRONGLY advised to place patients in a safe environment and to stay close to them throughout the session.

It is also recommended to increase the duration and intensity of the stimulation very gradually, after an initial short session to make sure of patients' tolerance to this type of stimulation.

Virtualis declines any liability for any disorders suffered by patients during or after use of its software.

1.2. Hardware and minimum configuration requirements

Hardware required to use the system:

- VR Ready PC
- VR System: HTC VIVE, HTC VIVE Pro or compatible system
- Lighthouse bases (HTC VIVE tracking)
- HTC VIVE Controller
- XBOX 360 Controllers
- USB HUB

In order to install and use our virtual reality applications, we recommend a configuration equal to or higher than the system requirements:

Technical Minimum Requirements GPU NVIDIA: Gen9 GTX 970 / Gen10 GTX 1060 AMD Radeon: R9 290 / RW 480 / Vega 56 Operating System Windows 7 SP1 RAM 8 Go

VIRTUALIS Page **3** of **10**




2. USE OF PATIENT MANAGEMENT

Once connected to the Patient Management software, you arrive on the home page. It is from this home page that you will be able to start your VR software as well as the other Patient Management functions.

The software can be grouped according to criteria such as "Assessment" or "Re-education" and then by pathology type: Neurology, Balance, Functional or Travel sickness.


You can start or switch from one software to another from the home page by clicking the corresponding "Start" or "Protocols" button.




Some software can be started either in *manual mode*, by directly clicking the "Start" button, or in *protocol mode* by clicking the "Protocols" button.

The *manual mode* allows users to choose the settings for each environment. The *protocol mode* offers several sessions with different difficulty levels to test and gradually accustom patients to the VR environment.

VIRTUALIS Page 4 of 10



Software that is not part of your subscription package is greyed out. If you want to use it, please contact our sales department.



VIRTUALIS Page **5** of **10**



3. MEMORIZATION

3.1. Start interface




When starting the software in *manual mode* ("Start" button) from Patient Management, the software is opened in a start interface consisting of a module selection menu on the left, a set up area on the right, and an action area at the bottom right.

Depending on the module selected in the left menu, the set up area shows the various possible settings/information.

The general Patient Management menu can be accessed from the start interface by simply clicking the "Back" button located in the action area, or by pressing the "escape" key on the keyboard.

The module is launched by simply clicking the "start" button in the action area.



Once this button is pressed, the module starts by taking into account the specified settings.

The Start/Quit buttons allow the environment to be played back or stopped entirely to adapt the experience to the patient's sensations.

VIRTUALIS Page 6 of 10



3.2. Module field of application

Memorization exercises for cognitive impairments

3.3. Installing the patient

Patients can be seated for the "Memorization" module.

3.4. Session settings

The variable settings for this module are as follows:


Timer

Time given to patients to memorize the objects to be found.



Number of visible products

Uses the cursor to select the number of products displayed on the screen at the same time. The fewer products there are, the easier it will be for users to find the targeted products.




Selected articles

Selection of the items that patients will have to find in order. To select an object, drag it to the left area or double-click it. Random object selection can be used using the "Random" button by specifying the number of products to be selected.

VIRTUALIS Page **7** of **10**





Patient interface

Users can use the control to select the required articles and then select the 'Finish' button to confirm their choices.



VIRTUALIS Page 8 of 10



Score

At the end of the exercise, users will get a score showing their success relative to the selected articles.




3.5. Shortcuts

Keyboard, joystick and controller shortcuts can be accessed in two ways:

- on the "Shortcuts" tab available at the start interface level
- within the module, by clicking on the joystick icon in the upper right corner of the screen

VIRTUALIS Page **9** of **10**







3.6. Data processing

Data retrieval and analysis uses the Patient Management software.

VIRTUALIS Page 10 of 10