

SPATIAL NAVIGATION REF



Class I Medical Device

User manual

Distribution mode

Available for direct download at http://virtualisvr.com/espace-client/ Use under licence



Avenue de l'Europe - 34830 CLAPIERS - Tel. 09 80 80 92 91



DESCRIPTION

SPATIAL NAVIGATION software is an immersive 3D simulation based on virtual reality technology, i.e. it allows a person to be immersed in an artificial digitally created world. **Spatial Navigation** is a software used for reintegrating Activities of Daily Living, memory training and executive function training (movement strategy, planning).

INDICATIONS

Assessment and rehabilitation of cognitive disorders (memory, executive functions).

CONTRAINDICATIONS

Epileptic patients, children under 15 years of age, pregnant women.

FOR USE BY

Healthcare professionals: Physiotherapists; Occupational therapists; Neuropsychologists; ENT doctors; Neurologists; PMR doctors (physical medicine and rehabilitation), etc.

Research Centers: CNRS, CHU, INSERM, etc.

WARNINGS AND CAUTIONS

During sessions, stay close to the patient in order to anticipate any loss of balance or discomfort caused by the use of virtual reality.

Define a working area of about $3m^2$ to allow for risk-free movements.

Take a 10 to 15 minute break every 30 minutes of use.

Potential adverse effects are those due to software, i.e. vomiting, malaise, dizziness, syncope

The accessories required to use the software may emit radio waves that can interfere with the operation of nearby electronic devices. If you have a pacemaker or other implanted medical device, do not use the product until you have taken advice from your doctor or the manufacturer of your medical device.





Table of Contents

1.	GEN	IERAL	3
	1.1.	Advice for use	3
	1.2.	Hardware and minimum configuration requirements	3
2.	USE	OF PATIENT MANAGEMENT	4
3.	SPA	TIAL NAVIGATION	6
	3.1.	Start interface	6
	3.2.	Module field of application	7
	3.3.	Installing the patient	7
	3.4.	Session settings	7
	3.5.	Shortcuts	9
	3.6.	Data processing1	10



1. GENERAL

1.1. Advice for use

This type of re-education must be undertaken progressively, especially in Virtual Reality where the stimulation is much more "powerful" than the traditional optokinetic stimulators.

These stimulations have the potential to cause certain disorders: Vasovagal syncope, epileptic seizures, migraines, etc. (Despite a test phase on more than 2000 patients. Similarly to previous generation optokinetics, caution is required)

The contraindications are identical: Mainly epilepsy and migraines.

As postural reactions can be spectacular, it is VERY STRONGLY advised to place patients in a safe environment and to stay close to them throughout the session.

It is also recommended to increase the duration and intensity of the stimulation very gradually, after an initial short session to make sure of patients' tolerance to this type of stimulation.


Virtualis declines any liability for any disorders suffered by patients during or after use of its software.

1.2. Hardware and minimum configuration requirements

Hardware required to use the system:

- VR Ready PC
- VR System: HTC VIVE, HTC VIVE Pro or compatible system
- Lighthouse bases (HTC VIVE tracking)
- HTC VIVE Controller
- XBOX 360 Controllers
- USB HUB

In order to install and use our virtual reality applications, we recommend a configuration equal to or higher than the system requirements:






2. USE OF PATIENT MANAGEMENT

Once connected to the Patient Management software, you arrive on the home page. It is from this home page that you will be able to start your VR software as well as the other Patient Management functions.

The software can be grouped according to criteria such as "Assessment" or "Re-education" and then by pathology type: Neurology, Balance, Functional or Travel sickness.


You can start or switch from one software to another from the home page by clicking the corresponding "Start" or "Protocols" button.




Some software can be started either in *manual mode*, by directly clicking the "Start" button, or in *protocol mode* by clicking the "Protocols" button.

The *manual mode* allows users to choose the settings for each environment. The *protocol mode* offers several sessions with different difficulty levels to test and gradually accustom patients to the VR environment.






Software that is not part of your subscription package is greyed out. If you want to use it, please contact our sales department.





3. SPATIAL NAVIGATION

3.1. Start interface



When starting the software in *manual mode* ("Start" button) from Patient Management, the software is opened in a start interface consisting of a module selection menu on the left, a set up area on the right, and an action area at the bottom right.

Depending on the module selected in the left menu, the set up area shows the various possible settings/information.

The general Patient Management menu can be accessed from the start interface by simply clicking the "Back" button located in the action area, or by pressing the "escape" key on the

keyboard.

The module is launched by simply clicking the "Start" button in the action area.

		Set	-up area	Version 1.5.6
Setting the session	Mode	Space navigation	~	
Manual				
Shortcuts	Timer 20 second	•		
Settings	Create a character every : 1	seconds	Ray of the dodge area of the user :	0.2 meters
Module selection				
menu	Неф		Proview	
	Bet the timer. Beliect the colors of your choice. The patient : - has X seconds (timer) to memorize the s - mult gin the color sections in the correc The practitioner validates when the patient h	as completed.		51 22
	Your VR operating system :	or	n area	Start 🕐 Quit 🌩

Once this button is pressed, the module starts by taking into account the specified settings. You also have the possibility to modify some settings when the module has been launched, using the mouse.

The Start/Quit buttons allow the environment to be played back or stopped entirely to adapt the experience to the patient's sensations.



3.2. Module field of application

Motor coordination, spatial navigation, decision-making

3.3. Installing the patient

Patient standing if possible.

3.4. Session settings

The variable settings for this module are as follows:

Timer

Time given to patients to memorize the colors to find.

Timer	60	seconds		

Character creation


Used to select the frequency at which characters appear

Dodge area radius

Used to determine the distance between the user and the other characters

List of aisles to browse


Choice of the aisles that patients will have to go through in order. To select an aisle color, drag it to the left area or double-click it.






Patient interface

Users can use the controller in their environment to reach the right color aisles.



Score

At the end of the exercise, users will get a score showing their success relative to the selected articles.





3.5. Shortcuts

Keyboard, joystick or controller shortcuts can be accessed in two ways:

- on the "Shortcuts" tab available at the start interface level
- within the module, by clicking on the joystick icon in the upper right corner of the screen

Quit		
	Show / Hide interface	
	Show / Hide user li	ist
	Show / Hide bas	sket list avoidance area
Ecnes F1 F2 F3 F4 F5 F6	F7 F8 F9 F10 F11 F12	2 Inde Ard Paule Syst Dir
$\begin{bmatrix} n & 1 & 2 & 3 & 4 & 5 & 6 & 7 \\ 2 & \dot{e} & & \# & , & f & f & -1 & \dot{e} \\ \end{bmatrix}$	8 9 ↓ 1 € + ↓ 1 + ↓ + ↓ + ↓ + ↓ + ↓ + ↓ + ↓ + ↓ +	Inser R 1
		Suppr Fin ▼ 7 8 9 L T L T L 1 L
	JKLM%	
];;;;;;} ;;∴;;;	
	Altigr 🖉 🚺 Ctri	
Recenter view	View Score finish	Increase characters
Show / Hide FPS	chara	number acters Jumber
		Decrease user avoidance area







Fast teleportation: Use the controller to aim at a sign to teleport to the associated aisle. Touch an item using the controller and then press the button to select it

3.6. Data processing

Data retrieval and analysis uses the Patient Management software.